
SPIRO[®]

SPEAKING VALVE

REF 710 CE 0402

REF 710R CE 0402

REF 711 CE

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en Instructions for Use

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


no Bruksanvisning

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sv Bruksanvisning

Symbol Glossary

Symbol	Title	Description
	Caution!	Indicates the need for the user to consult the Instructions for Use for important cautionary information, such as warnings and precautions.
	Consult instructions for use	Indicates the need for the user to consult the Instructions for Use.
	Manufacturer	Indicates the medical device manufacturer.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Use-by date	Indicates the date after which the medical device is not to be used.
	Batch code	Indicates the manufacturer's batch code so the batch or lot can be identified
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
	Keep dry	Indicates a medical device that needs to be protected from moisture.

Symbol	Title	Description
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use	Indicates a medical device that is intended for one use. (Article 740 only)
	Quantity	Indicates the total number of devices in the package.
QTY:	Quantity	Indicates the total number of devices in the package.
Rx Only	Prescription only (USA)	Indicates that the product is a medical device as defined in 21 CFR 820.3(l) and Federal Law (USA) restricts this device to sale by or on the order of a physician (21 CFR 801.109)
	CE mark	Indicates the manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

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SPIRO SPEAKING VALVE **REF** 710 + 710R + 711

Single patient use only

The SPEAKING VALVE is a one-way valve which opens ONLY FOR INHALATION. The exhalation passages MUST BE OPEN. Measurement of the tracheal pressure at passive exhalation with the Speaking Valve in place is advised (< 5 cm H₂O is recommended).

N.B. THE SPEAKING VALVE MUST NEVER BE USED IN COMBINATION WITH AN INFLATED CUFFED CANNULA WITHOUT FENESTRATION. ENSURE THAT THE PATIENT CAN TOLERATE A CANNULA WITH DEFLATED CUFF.

CONTRAINDICATIONS: This device is not for use with patients who have respiratory infection, extensive phlegm, stenosis, vocal cord paralysis or acute airway obstructions. This device is not for use with inline ventilators, endotracheal tubes or other artificial airways.

The COVER (5) is provided with a tapered section which is seated in the valve HOUSING (3) during use. The further pressed in, the tighter the cover is secured. Avoid pressing the cover too hard, as it reduces the valve's life and makes it difficult to remove.

DAILY CHECKS: Check valve operation daily and inspect for damage. Replace valve with new unit in the case of damage or defective operation. NEVER DISASSEMBLE THE MEMBRANE!

CLEAN the valve daily. Do not use sharp tools. Use the fingers or PLASTIC TWEEZERS (P/N 705, sold separately) to remove/re-insert the filter. The valve housing and cover shall be air dried after washing with luke-warm water and mild fragrance-free detergent. Rinse the valve cover in both directions to remove all cleaning liquids. Other detergents or solvents must not be used. Always check the operation again after cleaning and check that the air passages are free. Do not wash and re-use the filter!

The valve should always be used with FILTER (4) (P/N 740, sold separately). THE FILTER SHOULD BE REPLACED daily or more frequently if required. Ensure that the filter is seated evenly against the bottom of the valve housing, and that air flow through the filter is not obstructed by phlegm or other substances.

A HEAT SHIELD/SHOWER PROTECTOR (6) (P/N 730, sold separately) should be used as required with the inhalation opening facing downward.

An ADAPTER (2) (P/N 707, 708, sold separately) is only required for connection to a cannula (1) without an ISO 15 connector (taper).

SUPPLEMENTAL OXYGEN can connect via the Ø5.5 mm nipple on the valve housing. The WHITE O2 NIPPLE CAP (7) (P/N 710/710R only) shall be fitted to the Speaking Valve and removed only when oxygen is being supplied.

The SPEAKING VALVE should be replaced with a new unit after 29 days of use, or more frequently if necessary. The valve must not be used while the patient is sleeping or otherwise when the patient's consciousness is impaired. Use a SPIRO HUMIDIFIER (P/N 702, sold separately) in such cases.

The Speaking Valve must be removed before inhalation of medicine.

RESISTANCE TO FLOW for SPIRO Speaking Valve is < 2 cm H₂O*, and with the heat shield/shower protection on < 3 cm H₂O*.

All products should be stored in a dark and dry environment and at normal indoor temperature.

For the STORAGE PERIOD OF SPIRO PRODUCTS: See year and month next to the 'HOURGLASS' symbol on the product label. The products shall be disposed of in accordance with the applicable national regulations.

ALWAYS CONSULT YOUR DOCTOR, SUPPLIER OR MANUFACTURER IF YOU HAVE ANY QUESTIONS.

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Fogless International AB