


SPIRO[®] Catalogue 2015







*Improving the quality of life
for tracheostomy patients*

F  **GLESS**
International

HME's	2
SPEAKING VALVES	3
BASIC SETS	5
ACCESSORIES	6
ORDERING NOTES	7





HME's

	REF	DESCRIPTION	QTY
	702	SPIRO 702 Humidifier is a HME for spontaneous breathing tracheostomy patients. This device provides comfortable heat and moisture to the patient and prevents coarse particles from entering the respiratory system. The device has an exchangeable filter. The main device housing can be re-used multiple times with correct care. The patient cannot speak with this device and is it used mostly when the patient is resting or sleeping.	
	702-2	2 x SPIRO 702 Humidifier packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Bag dimensions: 118 x 70 mm • Total package weight: 15 gr. • Minimum Order: 10 devices • 1 x Instructions For Use per 10 devices 	2
	702-2K	2 x SPIRO 702 Humidifiers individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 83 x 72 x 42 mm • Total package weight: 35 gr. • Minimum Order: 20 devices (10 cartons) • 1 x Instructions For Use per carton 	2
	702-5K	5 x SPIRO 702 Humidifiers individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 157 x 72 X 42 mm • Total package weight: 67 gr. • Minimum Order: 50 devices (10 cartons) • 1 x Instructions For Use per carton 	5



SPEAKING VALVES

	REF	DESCRIPTION	QTY
	710	SPIRO 710 Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients. Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent coarse particles entering the respiratory system. The main valve housing can be re-used multiple times with correct care. When resting or sleeping the patient should use the SPIRO 702 Humidifier. This transparent device is recommended for use in the home care environment.	
	710-1	1 x SPIRO 710 Speaking Valve packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Bag dimensions: 118 x 70 mm • Total package weight: 6 gr • Minimum Order: 10 devices • 1 x Instructions For Use per 10 devices 	1
	710-2K	2 x SPIRO 710 Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 83 x 72 x 42 mm • Total package weight: 31 gr. • Minimum Order: 20 devices (10 cartons) • 1 x Instructions For Use per carton 	2
	710-5K	5 x SPIRO 710 Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 157 x 72 X 42 mm • Total package weight: 57 gr. • Minimum Order: 50 devices (10 cartons) • 1 x Instructions For Use per carton 	5
	710R	SPIRO 710R Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients. Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent coarse particles entering the respiratory system. The main valve housing can be re-used multiple times with correct care. When resting or sleeping the patient should use the SPIRO 702 Humidifier. This coloured device is recommended for use in the emergency or intensive care hospital environment.	
	710R-1	1 x SPIRO 710R Speaking Valve packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Bag dimensions: 118 x 70 mm • Total package weight: 6 gr • Minimum Order: 10 devices • 1 x Instructions For Use per 10 devices 	1






SPEAKING VALVES

	REF	DESCRIPTION	QTY
	710R-2K	2 x SPIRO 710R Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 83 x 72 x 42 mm • Total package weight: 31 gr. • Minimum Order: 20 devices (10 cartons) • 1 x Instructions For Use per carton 	2
	710R-5K	5 x SPIRO 710R Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 157 x 72 X 42 mm • Total package weight: 57 gr. • Minimum Order: 50 devices (10 cartons) • 1 x Instructions For Use per carton 	5
	711	SPIRO 711 Speaking Valve without oxygen port for spontaneous breathing tracheostomy patients. Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent coarse particles entering the respiratory system. The main valve housing can be re-used multiple times with correct care. When resting or sleeping the patient should use the SPIRO 702 Humidifier. This transparent device is recommended for use in the home care environment.	
	711-1	1 x SPIRO 711 Speaking Valve packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Bag dimensions: 118 x 70 mm • Total package weight: 6 gr • Minimum Order: 10 devices • 1 x Instructions For Use per 10 devices 	1
	711-2K	2 x SPIRO 711 Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 83 x 72 x 42 mm • Total package weight: 31 gr. • Minimum Order: 20 devices (10 cartons) • 1 x Instructions For Use per carton 	2
	711-5K	5 x SPIRO 711 Speaking Valves individually packaged in a 300 gsm eco-friendly paper carton. <ul style="list-style-type: none"> • Dimensions: 157 x 72 X 42 mm • Total package weight: 57 gr. • Minimum Order: 50 devices (10 cartons) • 1 x Instructions For Use per carton 	5

BASIC SETS

	REF	DESCRIPTION	QTY
Under Development	TBA	<p>SPIRO Complete Set Speaking Valve – Including 1 x SPIRO 710R Speaking Valve, 1 x 702 Humidifier, 1 x 740-25 Filter Sheet, 1 x 705 Tweezers, 1 x 703 Shower Cap, 1 x 730 Shower Cap.</p> <ul style="list-style-type: none"> • Packaging: PET • Dimensions: • Total package weight: • Minimum Order: 10 sets • 1 x Instruction For Use per set 	1
	760	<p>SPIRO Basic Set Speaking Valve – Including 1 x SPIRO 710 Speaking Valve, 1 SPIRO Humidifier, 1 x 705 Tweezers.</p> <ul style="list-style-type: none"> • Packaging: PET • Dimensions: 150 x 130 x 45 mm • Total package weight: 56 gr. • Minimum Order: 10 sets • 1 x Instruction For Use per set 	1
	770	<p>SPIRO Basic Set Humidifier – Including 2 x 702 Humidifiers, 1 x 705 Tweezers,</p> <ul style="list-style-type: none"> • Dimensions: 150 x 130 x 45 mm • Total package weight: 56 gr. • Minimum Order: 10 sets • 1 x Instruction For Use per set 	1

ACCESSORIES

	REF	DESCRIPTION	QTY
	703	<p>SPIRO 703 Shower Cap for use with the SPIRO 702 Humidifier. Intended to prevent water entering the tracheostomy opening when the patient is showering. Can also be used to direct warm air from the body in cold climates. Transparent in colour.</p> <ul style="list-style-type: none"> • Packaging: 2 pieces per polyethylene zip lock bag • Minimum Order: 10 pieces 	1
	705	<p>Tweezers. For exchanging of the filter.</p> <ul style="list-style-type: none"> • Packaging: 2 pieces per polyethylene zip lock bag • Length: 130 mm • Minimum Order: 10 pieces 	1
	707	<p>Adaptor for silver cannulae 13/11 mm. Transparent in colour.</p> <ul style="list-style-type: none"> • Packaging: 2 pieces per polyethylene zip lock bag • Minimum Order: 10 pieces 	1
	708	<p>Adaptor for silver cannulae 14/11 mm. White in colour.</p> <ul style="list-style-type: none"> • Packaging: 2 pieces per polyethylene zip lock bag • Minimum Order: 10 pieces 	1
	730	<p>SPIRO 730 Shower Cap for use with the SPIRO 710, 710R and 711 Speaking Valves. Intended to prevent water entering the tracheostomy opening when the patient is showering. Can also be used to direct warm air from the body in cold climates. Transparent in colour.</p> <ul style="list-style-type: none"> • Packaging: 2 pieces per polyethylene zip lock bag • Minimum Order: 10 pieces 	1

ACCESSORIES

	REF	DESCRIPTION	QTY
	740	SPIRO Filter. Exchangeable filter which is compatible with the SPIRO 702 Humidifier and SPIRO 710, 710R and 711 Speaking Valves. The filter prevents coarse particles from entering respiratory system.	
	740-25	SPIRO Filter – 1 sheet containing 25 filters. Packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Dimensions: 170 x 170 x 10 mm • Minimum Order: 10 packs (250 filters) 	25
	740-100	SPIRO Filter – 4 sheets containing 25 filters each. Packaged in a polyethylene zip-lock bag. <ul style="list-style-type: none"> • Dimensions: 170 x 170 x 35 mm • Minimum Order: 10 packs (1000 filters) 	100

ORDERING NOTES

1. Stated order quantities should be on a per device or per piece basis.
2. All products are delivered ex-works from Vaggeryd, Sweden.
3. Delivery time 2-4 weeks after confirmation of order.
4. Terms of payment: Cash on delivery (letter of credit) or according to agreement.
5. Service fee to be added for final packaging and handling.
6. Fogless® is SS-EN ISO 9001:2008, SS-EN ISO 13485:2012 certified and all products are CE-marked (1 or 2a) in accordance with MDD 93/42/EEC. 710/711 FDA 510(k) approved.

SPIRO® Speaking Valve advantages:

- *low inhalation resistance*
- *filter to prevent coarse particles entering the respiratory system*
- *three different variations to suit the individual user*
- *unique accessories*
- *quiet operation and easy to maintain*
- *integrated oxygen nipple (710/710R)*
- *Swedish quality design and manufacture*

